

SmartMan QCI Module – 5.0 Export Data

v2.0 and later

Contents

SmartMan QCI Module – 5.0 Export Data	1
5.0 UNDERSTANDING EXPORTED DATA	2
5.1 KEY TO ACTIVITIES	2
5.2 DATA RECORDS	2
5.2.1 Individual Record	2
5.2.2 Category (Variable)	2
5.3 GROUP ONE HEADINGS - IDENTIFICATION	2
5.3.1 Key Code	2
5.3.2 Student Number	2
5.3.3 Surname	2
5.3.4 Activity	3
5.3.5 Date	3
5.4 GROUP TWO HEADINGS - SCORES	3
5.4.1 Score 1 = Compression Depths (C/CPR) or Inspiration Volume (V)	3
5.4.2 Score 2 = Compression Durations (C/CPR) or Inspiration Time (V)	3
5.4.3 Score 3 = Compressions Release (C/CPR) or Interval Score (V)	4
5.4.4 Score 4 = Average Score 1;2;3	4
5.4.5 Score 5 = Compliant Compressions [depth; rate; release] (C/CPR) or Compliant ventilations [volume; duration] (V) [weighted against ideal totals]	4
5.4.6 Score 6 = Cycle Time Score	4
5.4.7 Score 7 = Final Score (weighted Score 4;5;6)	4
5.4.8 Score 8 = Total Time/Vents in activity	5
5.4.9 Score 9 = Hands Off Time (Total) (C/CPR) or 0 (V)	5
5.4.10 Score 10 = Hands Off Time (exclude time for ventilations) (C/CPR) or 0 (V)	5
5.4.11 Score 11 = Duty Cycle (C/CPR) or 0 (V)	5
5.4.12 Score 12 = Compressions performed (C/CPR) or Ventilations performed (V)	6
5.4.13 Score 13 = Inspiration Volume (CPR) or 0 (C/V)	6
5.4.14 DD = Depth and Duration	6
5.4.15 DDR = Depth, Duration and Rate	6
5.5 GROUP THREE HEADINGS – TOTALS & OTHER	6
5.5.1 TotalC = Total compressions (C/CPR) or 0 (V)	6
5.5.2 TotalV = Total ventilations (CPR/V) or 0 (C)	6

5.5.3 Level = With or With Access to Feedback	7
5.5.4 File = File name of the results.	7
5.5.5 Professional Group = Assigned Group in QCI Module.....	7

5.0 UNDERSTANDING EXPORTED DATA

When you have filtered the results in the SmartMan QCI Module, and you are satisfied with what they are displaying, you may wish to export those filtered results.

The data in the exported file contains more detailed information. The exported data meets the filtered parameters that you have set in the QCI Module.

In Excel each record will contain a lot of information. Each column will have a heading. Below is a list of these headings and a detailed explanation of how the data in that column is derived.

5.1 KEY TO ACTIVITIES

Each CPR activity on SmartMan can produce results.

C = Compression activities;
V = Ventilation activities)
CPR = CPR activities

5.2 DATA RECORDS

To understand what is contained in each category, column, see below. The bold faced words below will be the headings in each column. Although below the categories move down, in the spreadsheet the column headings move from left to the far column right.

5.2.1 Individual Record

The data results for the activity for each individual will be a single row in the Excel spread sheet.

5.2.2 Category (Variable)

Each column in that row will contain pre-processed data related to the heading of that column.

5.3 GROUP ONE HEADINGS - IDENTIFICATION

The first group of columns contain general identifying information.

5.3.1 Key Code

This is the key code issued by SmartMan. All results are related to a specific key code. This is the key code under which this set of results was performed. Different key codes can be used for different groups so that it is easy to identify them. You can get more key codes from SmartMan.

5.3.2 Student Number

When a student is enrolled, the key code is linked to a unique identified. This unique identifier is the student number.

5.3.3 Surname

This is the last name of the person enrolled and under whose name these results were performed.

5.3.4 Activity

This is the practice or test activity performed on SmartMan. This identifies which activity was performed on SmartMan and thus which results you are looking at. Here is a list of the activities.

Abbreviations for Activities

CP/T = Compression Practice/Test

VPCPR/VTCP = Ventilations in CPR Practice/Test

RBP/T = Rescue Breathing Practice/Test

CPR1RP/T = CPR 1 Rescuer Practice/Test

CPR2RP/T = 2 Rescuer CPR Practice/Test (lay rescuer BMV mode)

CPRAAP/T = 2 Rescuer Advanced Airway Practice/Test

CPRATP/T = 2 Rescuer Advanced Airway; Timed Practice/Test

5.3.5 Date

This is the date that was recorded when the file was created. It depends on the date setting on the computer when the results were performed. It is created when the practice or the test is performed. No results can be displayed unless this file is created and the date is put onto it.

5.4 GROUP TWO HEADINGS - SCORES

The second group of headings refer to scores. Below is an explanation of each score and how they are derived.

5.4.1 Score 1 = Compression Depths (C/CPR) or Inspiration Volume (V)

This is a Percentage score.

For compressions or CPR activities, Score 1 represents the number of compressions performed at the correct depth, expressed as a percentage of the total compressions performed in the activity. The current values are 5 - 6 cm

In ventilation activities Score 1 represents the number of ventilations performed at the correct volume of inspiration expressed as a percentage of the total number of ventilations performed in the activity. The current volume is 0.5 to 0.7 L.

5.4.2 Score 2 = Compression Durations (C/CPR) or Inspiration Time (V)

This is a Percentage score.

For compressions or CPR activities, Score 2 represents the number of compressions performed at the correct rate, expressed as a percentage of the total number of compressions performed in the activity. The current rate is 100 - 120 compressions per minute.

For ventilation activities, Score 2 represents the number of ventilations performed at the correct duration of inspiration, expressed as a percentage of the total number of ventilations performed in the activity. The current duration of inspiration is 1 s.

5.4.3 Score 3 = Compressions Release (C/CPR) or Interval Score (V)

This is a Percentage score.

For compressions or CPR activities, Score 3 represents the number of compressions fully released expressed as a percentage of the total number of compressions performed in the activity.

For ventilation activities, Score 3 represents the number of intervals between ventilations (that correspond to the correct rate of ventilations), expressed as a percentage of the total number of ventilations performed in the activity. For rescue breathing at 10 - 12 breaths per minute, intervals should be 5 - 6 seconds between the start of successive ventilations.

5.4.4 Score 4 = Average Score 1;2;3

This is derived by taking the total scores in score 1 + score 2 + score 3 and finding the average.

5.4.5 Score 5 = Compliant Compressions [depth; rate; release] (C/CPR) or Compliant ventilations [volume; duration] (V) [weighted against ideal totals]

This is a Percentage score.

For compressions or CPR activities, Score 5 represents the number of individual compressions performed at the correct rate, correct depth and fully released, expressed as a percentage of the total number of compressions performed in the activity.

For ventilation activities, Score 5 represents the number of individual ventilations performed at the correct duration and the correct volume of inspiration, expressed as a percentage of the total number of ventilations performed in the activity. The current duration of inspiration is 1 s.

5.4.6 Score 6 = Cycle Time Score

This is a score a percentage.

For compressions or CPR activities, Score 6 represents the difference in the time the learner took to complete the activity compared with the ideal time to complete the activity. The ideal time is based on the median of the compression rate range (100 to 120 compressions per minute) and includes or excludes ventilation cycles (as appropriate).

For ventilation activities (Rescue Breathing), Score 6 represents the difference in the time the learner took to complete the activity compared with the ideal time to complete the activity.

For ventilation activities (CPR Ventilations), Cycle Time is not evaluated an a nominal 100% is allocated to each score

5.4.7 Score 7 = Final Score (weighted Score 4;5;6)

This score is the SmartMan Percentage which is shown on screen when an activity is completed. This is a weighted score to reflect that some value can be done for a victim even though a skill is performed slightly outside of the target performance range on various skill parameters. For example, if all compressions were very close to deep enough but slightly below what is specified in the guidelines, and all were performed at the correct rate and with full release, then a compliance score would be zero.

The SmartMan score provides the user with a consistent reasonable indication of the value of their performance to a potential victim. In an educational sense it also provides a student with a way to compare one performance with another to gauge whether they have improved or not and whether they can maintain their skill level.

Here are the current weightings to derive that percentage score.

Compression Activities

Score 7 = Score 4 (30%) + Score 5 (50%) + Score 6 (20%)

Score 7 is factored by the variance in number of compressions performed from the ideal number.

CPR Activities

Score 7 =

Compression Components = Score 4 (30%) + Score 5 (30%) + Score 6 (10%) +

Ventilation Components = Score 4 (15%) + Score 5 (15%)

Score 7 is factored by the variance in number of compressions performed from the ideal number.

Ventilation Activities

Score 7 = Score 4 (30%) + Score 5 (40%) + Score 6 (15%) + Score 8 (15%)

5.4.8 Score 8 = Total Time/Vents in activity

In Compression activities, score 8 is the total time (s) to complete the whole activity compared to the ideal time.

In Rescue Breathing Activities, score 8 is a percentage of the number of ventilations actually performed compared to the ideal (10 to 12) ventilations. For example many people perform 25 to 35 rescue breaths per minute. This would show up in this score.

In the practice of Ventilations for CPR activity, it is expecting a minimum of 5 pairs of ventilations, above that scores 100%

5.4.9 Score 9 = Hands Off Time (Total) (C/CPR) or 0 (V)

This is the Total Time for when your hands are off the chest. This includes time between compressions and it includes the time off the chest to perform ventilations during a CPR activity.

This can be important if a person is taking a very long time to achieve a ventilation. The target time for CPR is to give 2 vents in 4 seconds. If it takes 5 or 6 seconds to get a ventilation, it will show up here.

5.4.10 Score 10 = Hands Off Time (exclude time for ventilations) (C/CPR) or 0 (V)

This is the time the Hands are off of the chest but does not include the time for ventilations. This number allows you to pick up poor performances in how the chest compressions are done. For example with a person jabbing at the chest, this number is too high. A chest compression is supposed to be a smooth continuous performance. Ventilation times are not recorded in this number.

5.4.11 Score 11 = Duty Cycle (C/CPR) or 0 (V)

Duty Cycle appears in the latest guidelines. Research has shown that not all compressions are equal in creating forward flow of blood through the brain. The research states that you are creating much better pressure to move the blood forward as you lower the duty cycle score.

In terms of how this relates to the actual performance, chest compression should be a smooth continuous down and release. The research suggests that a score below 50% is helpful. Excellent scores are 35% and lower

This is the amount of time you spent compressing (pushing down plus releasing) during a cycle.

“Duty Cycle is important as “Coronary blood flow is determined partly by the duty cycle.”
Guidelines 2010.

5.4.12 Score 12 = Compressions performed (C/CPR) or Ventilations performed (V)

This is a data dump of the cycle array.

For compressions and CPR activities, Score 12 represents the number of individual compressions performed in each cycle of the activity. Each cycle is delimited by a "|".

For ventilations activities, Score 12 represents the number

5.4.13 Score 13 = Inspiration Volume (CPR) or 0 (C/V)

This score is the percentage of ventilations that were performed at the correct volume range of 0.5L to 0.7L. In ventilation only activities this would be Score1, but in CPR activities Score1 is the percent of compressions performed at the correct depth, hence this extra column is required when reporting for ventilations only activities.

5.4.14 DD = Depth and Duration

This relates only to chest compressions. D=Depth of the compression performed and the second D=Duration of the compression performed (this is rate at which it is performed).

This score is the raw number of compressions that were performed in compliance with the current guidelines for D and D. That is for every compression where both Depth and Duration is performed within the guidelines, it is counted one.

This number is a tally of the number of compressions performed correctly for both D and D. Thus this number tells us the number of compressions that were compliant for both Depth and Duration.

5.4.15 DDR = Depth, Duration and Rate

This relates only to chest compressions. D=Depth of the compression performed and the second D=Duration of the compression performed (this is rate at which it is performed) and , R=Release of the chest for the compression.

This score is the raw number of compressions that were performed in compliance with the current guidelines for D and D and R. That is, for every compression where Depth, Duration and Rate are performed within the guidelines, it is counted one.

This number is a tally of the number of compressions performed correctly for D, D and R. Thus this number tells us the number of compressions that were fully compliant.

Thus if you take the DDR number and divide it into the expected number of compressions for that activity, you will get a score for the percentage score of how many compression were performed correctly, with no weighting according to the guidelines. This gives you the compliance percentage.

5.5 GROUP THREE HEADINGS – TOTALS & OTHER

The third group of headings will provide totals, what level the activity was performed at, and it will contain a file identifier for the profession if that has been created in the QCI Module.

5.5.1 TotalC = Total compressions (C/CPR) or 0 (V)

This is a count. It is the total number of compressions performed. It is common to have people not perform exactly 30 in each cycle. This gets picked up here.

5.5.2 TotatV = Total ventilations (CPR/V) or 0 (C)

This is a count. It is the total number of ventilations performed. It is common for people to perform numbers different from those specified in the guidelines for various activities

5.5.3 Level = With or With Access to Feedback

The level indicates whether the real time feedback was displayed on the screen or not. If the feedback is being displayed, even if the computer is facing away from the subject, it is possible for coaching to occur.

S = Standard. This means that the person performing the activity and/or a supervisor had access to viewing the colored bars as the skills were performed.

A = Advanced. This means that the colored bars were not displayed when the activity was performed. You can be certain that these results were done as a blind activity.

5.5.4 File = File name of the results.

This is the original file name that was created when the activity was performed. It contained the original raw data points.

5.5.5 Professional Group = Assigned Group in QCI Module

In the QCI Module you are able to assign specific results to a professional group such as paramedic, nurse, etc. This will appear in this column. It is a convenient way to sort groups to examine specific performance parameters for that group.