

INSTALL SMARTMAN SOFTWARE AND USB DRIVERS

1

- Do **NOT** Plug In USB Cable Yet
- **Download** Software (<http://www.ambusmartman.com/s-update.html>)
- **Install** Software: Double click file. Follow prompts on screen. It will create directories, install Software, extract USB drivers, create a shortcut on your desktop.

2

- Plug in USB cable
- Choose option "Browse for driver software". Do **NOT** allow windows to install automatically. Point to the driver files which will be located here:
- **C:\Program Files\Ambu SmartMan v5.3\SmartManDriver,** OR
- **C:\Program Files (x86)\SmartManv5.3\SmartManDriver**

3

- Click on the SmartMan Desktop Icon to Start Program